




## ACKERMANN EUROPCLIP®

Clip Applicators & Titanium Atraumatic Haemostatic Clips



EUROPCLIP®

Clip Appliers

## EUROPCLIP® CLIP APPLIERS



- ▶ precise loading of the clips
- ▶ firm grip of the clips
- ▶ safe closing
- ▶ color code
- ▶ full range of sizes



Large:

<b>AI-321550</b>	Clip Appliers, straight, 190 mm
<b>AI-321552</b>	Clip Appliers, right angled, 190 mm
<b>AI-321554</b>	Clip Appliers, straight, 280 mm
<b>AI-321556</b>	Clip Appliers, right angled, 280 mm



Medium Large:

<b>AI-321540</b>	Clip Appliers, straight, 190 mm
<b>AI-321544</b>	Clip Appliers, right angled, 190 mm
<b>AI-321542</b>	Clip Appliers, straight, 280 mm
<b>AI-321546</b>	Clip Appliers, right angled, 280 mm



Medium:

<b>AI-321524</b>	Clip Appliers, straight, 150 mm
<b>AI-321526</b>	Clip Appliers, straight, 190 mm
<b>AI-321520</b>	Clip Appliers, right angled, 190 mm
<b>AI-321522</b>	Clip Appliers, straight, 280 mm
<b>AI-321528</b>	Clip Appliers, right angled, 280 mm



Small:

<b>AI-321530</b>	Clip Appliers, straight, 150 mm
<b>AI-321532</b>	Clip Appliers, straight, 190 mm
<b>AI-321534</b>	Clip Appliers, right angled, 190 mm



# EUROPCLIP® TITANIUM ATRAUMATIC HAEMOSTATIC CLIPS



## easy and convenient

- ▶ soft loading system
- ▶ sticky tape at the back of the cartridge

## safe and up-to-date

- ▶ the chevron shape facilitates the safe „tips-first“ closing action
- ▶ serrated inner surface and a heart shaped cross-section of the clips allow secure closing
- ▶ solid titanium gives a scatter-free CT image and a MIR compatibility

## sterile clips

- ▶ sterile packed clips in a peel pack



Large:

<b>AI74CL6</b>	6	20	120
	Clips/Cartridge	Cartridge/Box	Clips/Box



Medium Large:

<b>AI74CML6</b>	6	20	120
	Clips/Cartridge	Cartridge/Box	Clips/Box



Medium:

<b>AI74CM6</b>	6	20	120
	Clips/Cartridge	Cartridge/Box	Clips/Box



Small:

<b>AI74CS6</b>	6	20	120
	Clips/Cartridge	Cartridge/Box	Clips/Box



ACKERMANN INSTRUMENTE GmbH, GERMANY  
Eisenbahnstraße 65-67  
78604 Rietheim-Weilheim

Phone: +49 (0) 74 61 / 9 66 17-0  
Fax: +49 (0) 74 61 / 9 66 17-70

info@ackermanninstrumente.de  
www.ackermanninstrumente.de